
ANGIOJET[®]
U L T R A

Thrombectomy Set Information for Use

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Information for Use

Table of Contents

1. Device Description	2
Figure 1: AngioJet Ultra AVX Thrombectomy Set	2
Table 1: Thrombectomy Set Information	2
2. Indications and Usage	2
3. Contraindications	2
4. Warnings and Precautions	2
5. Adverse Events	3
6. Clinical Studies	3
Table 2: Randomized Study Comparing AngioJet vs. Surgical Embolectomy in A-V Access Conduits	3
7. Patient Counseling Information	3
8. Clinician Use Information	3
8.1 Inspection Prior to Use – Precautions	3
8.2 Materials Required (not included in the package)	4
8.3 Catheter Preparation	4
8.4 Catheter Delivery and Operation	4
9. Patient Information	5
10. How Supplied	5

1. DEVICE DESCRIPTION

The AngioJet Ultra Thrombectomy Set (See Figure 1) is one component of the AngioJet Ultra Thrombectomy System (AngioJet System or System). The other component, packaged and sold separately, is the multiple-use AngioJet Ultra Console (Console). The AngioJet Ultra Thrombectomy Set can only be used in conjunction with the AngioJet Ultra Console.

The integral components of the disposable Thrombectomy Set are:

- Catheter
- Pump
- Saline delivery tubing and waste tubing
- Collection bag


Figure 1: AVX Thrombectomy Set

A bag of sterile saline (not included) supplies the pump with heparinized saline through the saline delivery tubing. The pump pressurizes the saline.

The Thrombectomy Set uses high-velocity saline to create a low-pressure zone at the catheter tip. Thrombus is drawn into the catheter where it is fragmented by the jets and evacuated from the body. The waste tubing transports the thrombus debris from the catheter to the collection bag for ultimate disposal.

The Thrombectomy Set is designed for removing thrombus from A-V access conduits. The catheter is designed to track over a 0.035" guide wire and through a guide catheter or sheath having an inner diameter of 0.086" or larger, with sufficient clearance to allow manual contrast injection, if desired. See Table 1: Thrombectomy Set Information.

Table 1: Thrombectomy Set Information

Product	AVX
Minimum vessel diameter:	3mm
Working length:	50cm
Maximum shaft diameter:	6F
Wire platform:	OTW
Wire compatibility:	0.035"
Guide compatibility: inner diameter:	8F (≥ 0.086")
Sheath compatibility:	6F
<i>Recommended maximum run times (minutes),</i>	
Total run time:	10
Run time with blood flow:	5

2. INDICATIONS AND USAGE

The AngioJet Ultra AVX Thrombectomy Set is intended for use with the AngioJet Ultra System in breaking apart and removing thrombus from A-V access conduits ≥ 3.0 mm in diameter.

3. CONTRAINDICATIONS

Do not use the Thrombectomy Set in patients:

- Who are contraindicated for endovascular procedures;
- In whom the lesion cannot be accessed with the guide wire.

4. WARNINGS AND PRECAUTIONS*

- The Thrombectomy Set has not been evaluated for treatment of pulmonary embolism. There are reports of serious adverse events, including death, associated with cases where the catheter was used in treatment of pulmonary embolism.
- The AngioJet Ultra System should be used only by operators who have received appropriate training in its installation and use.
- Use the Thrombectomy Set only with the multiple-use AngioJet Ultra Console.
- Possible complications associated with the use of the AngioJet Ultra System are similar to those associated with other endovascular procedures. A list of possible complications is included in Section 5, Adverse Events. Operation of the catheter may cause embolization of some thrombus and/or thrombotic particulate debris. Debris embolization may cause distal vessel occlusion, which may further result in hypoperfusion or tissue necrosis.
- Cardiac arrhythmias during catheter operation have been reported in a small number of patients. Cardiac rhythm should be monitored during catheter use and appropriate management, such as temporary pacing, be employed if needed.
- Do not use the AngioJet Ultra System in patients who have a non-healed injury due to recent mechanical intervention in the vessel to be treated to avoid further injury, dissection, or hemorrhage.
- Do not use the catheter in vessels smaller than 3.0mm in diameter which may increase risk of vessel injury.
- Obstructing lesions that are difficult to cross with the catheter to access thrombus may be balloon dilated with low pressure (≤ 2 atm). Failure to pre-dilate difficult-to-cross lesions prior to catheter operation may result in vessel injury.
- It is recommended that the catheter be withdrawn before contrast delivery using a power injector.
- Do not use if sterile package is opened or damaged. The Thrombectomy Set is intended for single use only. Do not resterilize, reprocess, or reuse.
- Visually inspect the Thrombectomy Set prior to use to ensure that no damage has occurred during shipment. Do not attempt to straighten or use the catheter if it is bent or kinked. Attempting to do so may result in catheter rupture. Do not use a damaged Thrombectomy Set for patient treatment.
- Systemic heparinization is advisable to avoid pericatheterization thrombus and acute rethrombosis. This is in addition to the heparin added to the saline supply bag.

Physician discretion with regard to the use of heparin is advised.

- Operation of the AngioJet System causes transient hemolysis which may manifest as hemoglobinuria. Table 1 lists maximum recommended run times in a flowing blood field and total operating time for each Thrombectomy Set. Evaluate the patient's risk tolerance for hemoglobinemia prior to the procedure. Consider hydration prior to, during, and after the procedure as appropriate to the patient's overall medical condition.
- Large thrombus burdens in peripheral veins may result in significant hemoglobinemia which should be monitored to manage possible renal or other adverse events.
- The catheter should be operated over a 0.035" guide wire. Attempting to use a larger guide wire will damage the catheter and the guide wire.
- Ensure that the catheter tip is submerged in heparinized saline during operation outside the body to avoid undesired wetting of the sterile field or operator contamination.
- Thrombotic debris/fluid flow exiting the catheter manifold via the waste tubing should be monitored continually during use. If the effluent fluid is primarily saline (clear in appearance), catheter operation should be stopped and the catheter should be repositioned. Clear fluid indicates that the catheter may be occlusive within the vessel or operating in an area cleared of clot. Prolonged operation under occlusive conditions may result in tissue injury.
- During the procedure, do not retract the guide wire into the catheter. If retraction of the guide wire into the catheter occurs, it may be necessary to remove both the catheter and guide wire from the patient in order to re-load the catheter over the guide wire.
- Do not pull the catheter against abnormal resistance. If increased resistance is felt when removing the catheter, remove the catheter together with the sheath or guide catheter as a unit to prevent possible tip separation.

*See also Section 8, Clinician Use Information.

5. ADVERSE EVENTS

Potential adverse events which may be associated with use of the AngioJet Ultra Thrombectomy System are similar to those associated with other interventional procedures and include, but are not limited to:

- abrupt closure of treated vessel
- acute myocardial infarction
- acute renal failure
- bleeding from access site
- cerebrovascular accident
- death
- dissection
- embolization, proximal or distal
- hematoma
- hemolysis
- hemorrhage, requiring transfusion
- hypotension/hypertension
- infection at the access site
- pain
- perforation
- pseudoaneurysm
- reactions to contrast medium
- thrombosis/occlusion
- total occlusion of treated vessel
- vascular aneurysm
- vascular spasm

6. CLINICAL STUDIES

The results of a prospective, randomized U.S. IDE clinical trial comparing AngioJet vs. surgical embolectomy in A-V access conduits are summarized in Table 2.

Table 2: Randomized Study Comparing AngioJet vs Surgical Embolectomy in A-V Access Conduits (F105/LF140 Catheter Models)

	AngioJet	Surgical Embolectomy
Number of Grafts Treated	63	55
Acute Procedure Success (Angiographic Evaluation by Core Lab)	56/57 (98%)	37/38 (97%)
Functional Procedure Success (Ability to Dialyze through Graft within 72 Hours)	47/31 (77%)	45/55 (82%)
Grafts Requiring Same Session Secondary Thrombectomy Treatment	10/63 (16%)	2/55 (4%)**
Primary Patency Rate (Actuarial) at 30 Days	47%	53%
Primary Patency Rate (Linear) at 30 Days*	22/54 (41%)	27/52 (52%)
Secondary Patency Rate (Linear) at 30 Days	22/54 (41%)	27/52 (52%)
Grafts Re-Thrombosing and/or receiving Secondary Intervention within 30 Days	35/57 (61%)	26/53 (49%)
Primary Patency Rate (Actuarial) at 180 Days	13%	21%
Primary Patency Rate (Linear) at 180 Days*	4/41 (10%)	6/35 (17%)
Secondary Patency Rate (Linear) at 180 Days	18/41 (44%)	18/35 (51%)
Grafts Re-Thrombosing and/or receiving Secondary Intervention between 30 and 180 Days	14/44 (32%)	18/38 (47%)

* Does not include patients who experienced a final event (except for rethrombosis) prior to the reporting interval.

** Does not include same session surgical revision.

7. PATIENT COUNSELING INFORMATION

The physician should consider the following points in counseling the patient about the AngioJet Ultra System:

- Discuss the risk/benefit issues for this particular patient, both for AngioJet use and for other interventional treatments likely to be employed.
- Discuss post procedure home-care and rehabilitation guidelines.
- Discuss alterations to current lifestyle following treatment, both short and long term.

8. CLINICIAN USE INFORMATION

Refer to the Operator's Manual included with the AngioJet Ultra Console for additional AngioJet Ultra System operating instructions.

8.1 Inspection Prior to Use – Precautions

1. Verify "Use By" date printed on the product label has not expired. Do not use the Thrombectomy Set if "Use By" date has expired, as sterility or function may be compromised.

2. Inspect sealed sterile package before opening. If packaging integrity has been compromised, do not use, reprocess or resterilize the contents, and contact a Possis Medical, Inc. representative.

8.2 Materials Required (not included in the package)

- AngioJet Ultra Console
- Appropriate 0.035 inch guide wire
- Heparinized normal saline (HepNS), in bowl for wipe-downs and Thrombectomy Set priming (2000 units per liter (U/L) suggested)
- HepNS for injection, in bag to be hung on Console, and to be connected to the Thrombectomy Set (5000 U/L suggested)

8.3 Catheter Preparation

1. Heparinize a bag of sterile, room-temperature saline at a suggested rate of 5000 units per liter of saline and mix contents (1.0 liter is recommended, but is not included with the System). Hang the saline bag on the saline bag hook at the top of the Console.

Note: Physician discretion with regard to the use of heparin is advised.

2. Power up the AngioJet Ultra console, the console drawer will open.
3. Remove the catheter from the package using standard sterile technique. Inspect the catheter prior to use to ensure that no damage has occurred during shipment such as bends and kinks. Do not use or straighten and use a severely kinked product; such use may result in rupture and/or vessel injury. Do not use the catheter if any defects are noted, and contact a Possis Medical, Inc. representative.
4. Sterile Technician: Remove the catheter and sufficient tubing for ease of use from the sterile package and inspect for damage. Hand the rest of the Thrombectomy Set to the nonsterile technician for installation in the Console.

Note: If a clamp is used to keep the catheter in the sterile field, ensure that the clamp does not deform the tubing.

5. Nonsterile technician: Take the pump (not the piston head) from the sterile technician and insert the pump into the Console. Ensure that the waste bag tubing aligns with the roller pump.
6. Remove the cap from the Thrombectomy Set bag spike and insert the spike into the bag.
7. Push the button to close the Console drawer.
8. Place the foot pedal within easy access of the physician. Choose a location that will minimize accidental activation.
9. Prime the catheter by completely submerging the tip in heparinized saline and pressing the foot pedal.
10. Continue priming until the time display reaches zero seconds.
11. Confirm system set-up is successfully completed by removing foot from foot pedal. Status panel displays READY and green icon is illuminated.

8.4 Catheter Delivery and Operation

Note: AngioJet Ultra System thrombectomy treatments can be performed in conjunction with definitive treatments of the residual stenosis such as balloon angioplasty or stent placement.

1. The catheter is intended to be introduced through percutaneous access. A vascular sheath should be used to minimize trauma at the access site. Deliver the catheter to the thrombosed section. Ancillary guide catheters used during the procedure should have an inner diameter of at least 0.086" to allow sufficient room for passage of the catheter. A hemostasis valve should be attached to the guide catheter to prevent blood loss, and should have sufficient size to pass and seal around the catheter without restricting catheter movement. Guide wires or other standard devices and methods should be used to aid in placement of the sheath and/or guide catheter.

Caution: Obstructing lesions that are difficult to cross with the catheter to access thrombus may be balloon dilated with low pressure (≤ 2 atm). Failure to pre-dilate difficult-to-cross lesions prior to catheter operation may result in vessel injury.

2. After guide wire is advanced to the treatment site, the AVX Catheter will be back-loaded and advanced over the guide wire.

Tru-Seal® Hub Operation:

3. Before loading, open the hemostasis valve by pressing on the introducer on the AVX hub. (Fig. 2)
4. Load the AVX catheter over the 0.035" guide wire, and advance until the wire exits the Tru-Seal hub.
5. Close the hemostasis valve by grasping the introducer and pulling back, you will hear an audible click. (Fig. 2) No adjustment to the Tru-Seal hemostasis valve is necessary. If leaking occurs confirm the introducer piece is fully extended to the closed position. (Fig. 2)


Figure 2. Tru-Seal Operation

Caution: Do not exchange the guide wire. The guide wire should extend at least 3cm past the catheter tip at all times. During the procedure, do not retract the guide wire into the catheter. If retraction of the guide wire into the catheter occurs, remove both the catheter and guide wire from the patient in order to back-load the catheter over the guide wire. The catheter should always be back-loaded onto the guide wire to prevent the wire tip from exiting and binding in the catheter windows.

Caution: Do not pull the catheter against abnormal resistance. If increased resistance is felt when removing the catheter, remove the catheter together with the sheath or guide catheter as a unit to prevent possible tip separation.

6. Thrombectomy treatments may be accomplished by using the AngioJet Ultra System proximal to distal or distal to proximal. For peripheral A-V access conduit applications using the AVX Catheter, distal to proximal is recommended.

6.1 Catheter passes should be continuous and without deactivation of the AngioJet System. Repeated AngioJet System activation/deactivation is not advised, although several passes with the catheter may be made, as needed. Do not apply excessive force if advancement is difficult. **Ensure that the guide wire extends at least 3 cm past the catheter tip at all times. Do not retract the guide wire into the catheter during operation.**

6.2 Distal/proximal – Advance the guide wire across the thrombotic lesion. Advance the catheter over-the-wire to the distal end of the thrombotic lesion. Activate the AngioJet Ultra System by pressing the Console foot pedal and retract the catheter into the thrombus at a rate of approximately 3 mm/sec. At the proximal end of the thrombotic lesion, reverse direction to continue treatment in an antegrade manner.

6.3 Proximal/distal – Advance the guide wire across the thrombotic lesion. Advance the catheter over-the-wire, and position the catheter tip just proximal to the thrombus. Activate the AngioJet Ultra System by pressing the Console foot pedal and advance the catheter into the thrombus at a rate of approximately 3 mm/sec. At the distal end of the thrombotic lesion, reverse direction to continue treatment in a retrograde manner.

Note: The distal-to-proximal method should be used if the catheter can be positioned safely and when the thrombus is highly mobile with a proximal attachment point.

7. Standard contrast media, delivered by hand injection through the guide catheter or sheath and around the catheter, may be used to evaluate treatment. In some cases, it may be necessary to remove the catheter in order to inject a sufficient volume of contrast.

Caution: Operation of the AngioJet System causes transient hemolysis which may manifest as hemoglobinuria. Table 1 lists maximum recommended run times in a flowing blood field and total operating time for each Thrombectomy Set. Evaluate the patient's risk tolerance for hemoglobinemia prior to the procedure. Consider hydration prior to, during, and after the procedure as appropriate to the patient's overall medical condition.

Caution: Large thrombus burdens in peripheral veins may result in significant hemoglobinemia which should be monitored to manage possible renal or other adverse events.

Note: If the Thrombectomy Set is removed and/or inoperative, the Thrombectomy Set tubing, guide catheter and sheath should be flushed with HepNS to avoid thrombus formation and maintain lumen patency. Re-prime the Thrombectomy Set by submerging the tip in HepNS and operating it for at least 20 seconds before reintroduction to the patient.

8. Follow proper precautions for the handling of infectious waste. Reuse of Thrombectomy Set is prohibited. Refer to the Console Operators Manual for Console cleaning instructions and proper disposal of the Thrombectomy Set.

9. PATIENT INFORMATION

No separate patient information is provided because the catheter is used in association with other interventional procedures. Risks and benefits to the patient are similar to and, from the patient's perspective, a part of those interventional procedures.

10. HOW SUPPLIED

STERILE: The Thrombectomy Set is supplied sterile and non-pyrogenic. The Thrombectomy Set is sterilized with ethylene oxide gas (EO). Contents can be considered sterile in an unopened, undamaged package until the date shown on the product label.

CONTENTS: One (1) AngioJet Ultra Thrombectomy Set
One (1) *Information for Use* insert

EC
S
a
a
m
b
l
e

POSSIS®

Bringing Medical Possibilities to Life®

Manufactured and
Distributed by:

Possis Medical, Inc.
9055 Evergreen Boulevard N.W.
Minneapolis, MN 55433-8003
USA
tel +1.763.780.4555
toll free +1.800.810.7677
fax +1.763.780.2227

Authorized European
Representative:

CEpartner4U BV
Esdoornlaan 13,
3951 DB Maarn,
The Netherlands
tel +31 343 442 524
mobile +31 6 516 536 26
fax +31 343 442 162

AngioJet, Rheolytic, AVX, and Tru-Seal are trademarks of Possis Medical, Inc.,
Minneapolis, Minnesota

U.S. and foreign patents and/or patents pending

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